# **CRYOTUBES**

# CryoGen<sup>®</sup> Tubes / 1D linear barcode INTERNAL CAP Tubes





### CryoGen Tubes Internal cap

Sterile, DNase, RNase, Human DNA, ATP, Pyrogen free





Internal cap - case / 500

CAT. NO.	Vol. (ml)	Height (cm)	Bottom style	Internal packaging	Sale unit
CL1 ARBIPSTS	1.2	4.27	Self Standing	bag / 50	case / 500
CL2ARBIPSTS	2.0	4.82	Self Standing	bag / 50	case / 500
CL2ARIPSTS	2.0	4.70	Round bottom	bag / 50	case / 500
CL4ARBIPSTS	4.0	7.62	Self Standing	bag / 50	case / 500
CL4ARIPSTS	4.0	7.50	Round bottom	bag / 50	case / 500
CL5ARBIPSTS	5.0	9.03	Self Standing	bag / 50	case / 500
CL5ARIPSTS	5.0	8.92	Round bottom	bag / 50	case / 500

6 different colors caps available. To order add ending as follows: /AZ (blue) /BI(white) /G(yellow) /RO(red) /V(green) /VI(violet)





4.0 ml

Internal cap - packaging in bulk

CAT. NO.	Vol. (ml)	Height (cm)	Bottom style	Internal packaging	Sale unit
CL1ARBIPSTS/B	1.2	4.27	Self Standing	bag / 50	carton / 4500
CL2ARBIPSTS/B	2.0	4.82	Self Standing	bag / 50	carton / 4000
CL2ARIPSTS/B	2.0	4.70	Round bottom	bag / 50	carton / 4000
CL4ARBIPSTS/B	4.0	7.62	Self Standing	bag / 50	carton / 2500
CL4ARIPSTS/B	4.0	7.50	Round bottom	bag / 50	carton / 2500
CL5ARBIPSTS/B	5.0	9.03	Self Standing	bag / 50	carton / 2000
CL5ARIPSTS/B	5.0	8.92	Round bottom	bag / 50	carton / 2000



#### For automatic capper/decapper

Star shaped top cavity is engineered for use with automatic decapping equipment (see page 20)

## No contamination risk

Our innovative screw caps eliminate leakage and contamination. Both internal and external screw caps are co-molded with a thermoplastic elastomer (TPV) layer and are 95kPa certified to provide a 100% leak-proof seal. This proprietary molding process eliminates the risk of contamination associated with brands that use separate o-rings.

# Insert for color coding

Color coded inserts are available to fit both internal and external caps for quick sample identification (see page 24)





